

## EasyOne Pro Remote Service Firmware $\geq$ 1.3.4.0

### Procedure

Does your EasyOne Pro have unrestricted Internet access?


If yes, proceed with **Step A**

If no, proceed with **Step B**



### Step A Prepare for remote Service

1. Connected EasyOne Pro to the network
2. Switch on EasyOne Pro Utilities→Advanced→ Enter Password: EOPTM→Login→Remote Support→ Enter Password: EOPTM
3. Summit the ID-number and the Password to the ndd supporter



4. Troubleshoot
5. At the end: Close TeamViewer- Session by pressing Cancel