

Spirette™

Breathing mouthpiece for Infection Control



The ndd Spirette mouthpieces have been engineered to **fit perfectly into ndd devices to guarantee reliable and high quality measurement results.**

ndd Spirette mouthpieces **prioritise hygiene and safe re-use of ndd devices** by avoiding cross-contamination between patients.

All parts exposed to patient's breath are **single-use and disposable**

No expiration date

Technical Data

Material	Spirette: linear low density polyethylene (LLDPE) Bag: Polyethylene (PE)
Transport and storage information	Temperature: -25 to 70°C Rel. Humidity: 5 to 90%, non-condensing Atmospheric Pressure: 400 to 1060 hPa
Dimension (single Spirette)	25 x 38 x 140 mm (H x W x D), 8 g
Compatible with	EasyOne, Easy on-PC, EasyOne Pro and EasyOne Pro LAB

Order Information

Part number:	Product:
2050-1	Spirette, individually wrapped, dispenser box of 50 pcs.
2050-5	Spirette, individually wrapped, box of 200 pcs. (dispenser box included) wrapped, bag of 10 pcs.
2050-10	Spirette, individually wrapped, box of 500 pcs. (without dispenser box)