

# Spirette™

Breathing mouthpiece for Infection Control



The ndd Spirette mouthpieces have been engineered to **fit perfectly into ndd devices to guarantee reliable and high quality measurement results.**

ndd Spirette mouthpieces **prioritise hygiene and safe re-use of ndd devices** by avoiding cross-contamination between patients.

All parts exposed to patient's breath are **single-use and disposable**

No expiration date

## Technical Data

<b>Material</b>	Spirette: linear low density polyethylene (LLDPE) Bag: Polyethylene (PE)
<b>Transport and storage information</b>	Temperature: -25 to 70°C Rel. Humidity: 5 to 90%, non-condensing Atmospheric Pressure: 400 to 1060 hPa
<b>Dimension (single Spirette)</b>	25 x 38 x 140 mm (H x W x D), 8 g
<b>Compatible with</b>	EasyOne, Easy on-PC, EasyOne Pro and EasyOne Pro LAB

## Order Information

<b>Part number:</b>	<b>Product:</b>
2050-1	Spirette, individually wrapped, dispenser box of 50 pcs.
2050-5	Spirette, individually wrapped, box of 200 pcs. (dispenser box included) wrapped, bag of 10 pcs.
2050-10	Spirette, individually wrapped, box of 500 pcs. (without dispenser box)