



Medical Technologies

EasyOne Air

Portable & PC spirometer

Elevating spirometry. Improving outcomes.

n d d's **EasyOne Air** consistently leads the way in spirometry innovation. With a user-friendly touchscreen and integrated diagnostic prompts, this portable device seamlessly integrates into any practice workflow.

The advanced TrueFlow technology eliminates the necessity for calibration, ensuring that **EasyOne Air** produces readily actionable data to quickly inform decisions. When it comes to choosing a reliable spirometry solution, the clear and simple choice is **EasyOne Air**.

EasyOne



Simple. Certain. Proven.



Features

- Large color touchscreen
- Easy navigation
- Long-term stable calibration
- Rechargeable battery: >100 tests with one charge
- Bluetooth connectivity
- Flexible use
- Portable or PC mode
- Real-time animated incentives
- Real-time data transfer
- Proven integration with top EMR/EHR systems

Incentive screens

Available with *EasyOne Connect* software



Specifications

Test types: FVC, FVL, SVC, MVV and Pre/Post testing

Languages: Danish, Dutch, English, French, German, Italian, Norwegian, Polish, Portuguese, Russian, Spanish, Swedish

EMR connectivity: EasyOne Connect

Dimensions: 3.5" in diagonal size, 13 oz
88.9 mm in diagonal size, 0.356 kg

Standards: ATS/ERS 2022, 2019 and 2005, NIOSH, OSHA, SSA Disability

Article numbers

2500-2: EasyOne Air spirometry system

5050-50: EasyOne FlowTube box of 50

5050-200: EasyOne FlowTube box of 200

5050-500: EasyOne FlowTube box of 500
Not available in all countries



n d d Medical Technologies
300 Brickstone Square, Suite 604 Andover, MA, 01810 USA
Phone: +1 978 470 0923

n d d Medizintechnik AG
Technoparkstrasse 1, CH-8005 Zurich, Switzerland
Phone: +41 44 512 65 00



n d d's unique ultrasonic flow measurement is highly accurate in all flow ranges, independent of gas composition, pressure, temperature and humidity.

n d d **TrueFlow** is a resistance-free solution that does not require calibration during its lifetime.

EasyOne Connect

EasyOne Connect software offers one database and one platform for EasyOne point-of-care solutions. Data management has never been easier.



ID: 9001622 Rev: 04
© Copyright 2025. All rights reserved