



Medical Technologies

EasyOne Air

Portable & PC spirometer

Elevating spirometry.
Impacting patient
outcomes.

ndd's **EasyOne Air** consistently leads the way in spirometry innovation. With a user-friendly touchscreen and integrated diagnostic prompts, this portable device seamlessly integrates into any practice workflow.

The advanced TrueFlow technology eliminates the necessity for calibration, ensuring that **EasyOne Air** produces readily actionable data to quickly inform decisions. When it comes to choosing a reliable spirometry solution, the clear and simple choice is **EasyOne Air**.



EasyOne



Simple. Certain. Proven.



Incentive screens

Available with *EasyOne Connect* software



Features

- Large color touchscreen
- Easy navigation
- Long-term stable calibration
- Rechargeable battery: >100 tests with one charge
- Bluetooth connectivity
- Flexible use
- Portable or PC mode
- Real-time animated incentives
- Real-time data transfer
- Proven integration with top EMR/EHR systems

Specifications

Test types: FVC, FVL, SVC, MVV and Pre/Post testing

Languages: Danish, Dutch, English, French, German, Italian, Norwegian, Polish, Portuguese, Russian, Spanish, Swedish

EMR connectivity: EasyOne Connect

Dimensions: 3.4" x 6.1" x 1.4", 13 oz
87 mm X 155 mm x 36 mm, 0.356 kg

Standards: ATS/ERS 2022, 2019 and 2005, NIOSH, OSHA, SSA Disability

Article numbers

2500-2: EasyOne Air spirometry system

5050-50: EasyOne FlowTube box of 50

5050-200: EasyOne FlowTube box of 200

5050-500: EasyOne FlowTube box of 500
Not available in all countries

nnd
TrueFlow
makes the difference

nnd's unique ultrasonic flow measurement is highly accurate in all flow ranges, independent of gas composition, pressure, temperature and humidity.

nnd **TrueFlow** is a resistance-free solution that does not require calibration during its lifetime.

EasyOne Connect

EasyOne Connect software offers one database and one platform for EasyOne point-of-care solutions. Data management has never been easier.



nnd Medical Technologies
300 Brickstone Square, Suite 604 Andover, MA, 01810 USA
Phone: +1 978 470 0923

nnd Medizintechnik AG
Technoparkstrasse 1, CH-8005 Zurich, Switzerland
Phone: +41 44 512 65 00

