



Medical Technologies

# EasyOne Air

Portable & PC spirometer

Elevating spirometry.  
Impacting patient  
outcomes.

n d d's *EasyOne Air* consistently leads the way in spirometry innovation. With a user-friendly touchscreen and integrated diagnostic prompts, this portable device seamlessly integrates into any practice workflow.

The advanced TrueFlow technology eliminates the necessity for calibration, ensuring that *EasyOne Air* produces readily actionable data to quickly inform decisions. When it comes to choosing a reliable spirometry solution, the clear and simple choice is *EasyOne Air*.

EasyOne



# Simple. Certain. Proven.



## Incentive screens

Available with *EasyOne Connect* software



## Features

- Large color touchscreen
- Easy navigation
- Long-term stable calibration
- Rechargeable battery: >100 tests with one charge
- Bluetooth connectivity
- Flexible use
- Portable or PC mode
- Real-time animated incentives
- Real-time data transfer
- Proven integration with top EMR/EHR systems

## Specifications

**Test types:** FVC, FVL, SVC, MVV and Pre/Post testing

**Languages:** Danish, Dutch, English, French, German, Italian, Norwegian, Polish, Portuguese, Russian, Spanish, Swedish

**EMR connectivity:** EasyOne Connect

**Dimensions:** 3.4" x 6.1" x 1.4", 13 oz  
87 mm X 155 mm x 36 mm, 0.356 kg

**Standards:** ATS/ERS 2022, 2019 and 2005, NIOSH, OSHA, SSA Disability

## Article numbers

**2500-2:** EasyOne Air spirometry system

**5050-50:** EasyOne FlowTube box of 50

**5050-200:** EasyOne FlowTube box of 200

**5050-500:** EasyOne FlowTube box of 500  
Not available in all countries



n d d Medical Technologies  
300 Brickstone Square, Suite 604 Andover, MA, 01810 USA  
Phone: +1 978 470 0923

n d d Medizintechnik AG  
Technoparkstrasse 1, CH-8005 Zurich, Switzerland  
Phone: +41 44 512 65 00



**n d d**  
**TrueFlow**  
makes the difference

n d d's unique ultrasonic flow measurement is highly accurate in all flow ranges, independent of gas composition, pressure, temperature and humidity.

n d d **TrueFlow** is a resistance-free solution that does not require calibration during its lifetime.

## EasyOne Connect

**EasyOne Connect** software offers one database and one platform for EasyOne point-of-care solutions. Data management has never been easier.

SCAN ME



ID: LP-9001622 Rev: 05  
© Copyright 2025. All rights reserved