



Medical Technologies

Inspiring innovation. Every day.

EasyOne FlowTube

Breathing mouthpiece for
infection control



n d d **EasyOne FlowTube** mouthpieces avoid cross-contamination between patients and reduce cleaning of the device to an absolute minimum.

The **EasyOne FlowTube** mouthpieces have been engineered to fit perfectly into n d d's devices to guarantee reliable and high-quality measurement results.

All parts exposed to patient's breath are
single-patient-use and disposable.

Technical Data

Material

EasyOne FlowTube: polypropylene (PP)
Bag: low density polyethylene (LDPE)

Transport and storage information

Temperature: -25 to 70°C / -13° to 158°F
Rel. Humidity: 5 to 90%, non-condensing
Atmospheric Pressure: 400 to 1060 hPa

Dimensions (one EasyOne FlowTube)

21 x 38 x 150 mm, 11 g / 0.8 x 1.5 x 5.5", 0.38 oz

Compatible with

EasyOne Air, EasyOne Sky

Order Information

Part number:

Product:

5050-50

EasyOne FlowTube box of 50

5050-200

EasyOne FlowTube box of 200

5050-500

EasyOne FlowTube box of 500
Not available in all countries