



Medical Technologies

Inspiring innovation. Every day.

Spirette

Breathing mouthpiece for
infection control



n d d **Spirette** mouthpieces avoid cross-contamination between patients and reduce cleaning of the device to an absolute minimum.

The **Spirette** mouthpieces have been engineered to fit perfectly into n d d's devices to guarantee reliable and high-quality measurement results.

All parts exposed to patient's breath are
single-patient-use and disposable.

Technical Data

Material	Spirette: linear low density polyethylene (LLDPE) Bag: low density polyethylene (LDPE)
-----------------	---

Transport and storage information	Temperature: -25 to 50°C / -13 to 122°F Rel. Humidity: 5 to 90%, non-condensing Atmospheric Pressure: 400 to 1060 hPa
--	---

Dimensions (one Spirette)	25 x 38 x 140 mm, 10 g / 1 x 1.5 x 5.5", 0.35 oz
----------------------------------	--

Compatible with	Easy on-PC, EasyOne Pro, and EasyOne Pro LAB
------------------------	--

Order Information

Part number:	Product:
2050-1	Spirette box of 50
2050-5	Spirette box of 200
2050-10	Spirette box of 500 Not available in all countries