



Medical Technologies

Inspiring innovation. Every day.

Spirette

Breathing mouthpiece for infection control



n d d **Spirette** mouthpieces avoid cross-contamination between patients and reduce cleaning of the device to an absolute minimum.

The **Spirette** mouthpieces have been engineered to fit perfectly into n d d's devices to guarantee reliable and high-quality measurement results.

All parts exposed to patient's breath are **single-patient-use and disposable**.

Technical Data

| | |
|--|---|
| Material | Spirette: linear low density polyethylene (LLDPE) Bag: low density polyethylene (LDPE) |
| Transport and storage information | Temperature: -25 to 50°C / -13 to 122°F Rel. Humidity: 5 to 90%, non-condensing Atmospheric Pressure: 400 to 1060 hPa |
| Dimensions (one Spirette) | 25 x 38 x 140 mm, 10 g / 1 x 1.5 x 5.5", 0.35 oz |
| Compatible with | Easy on-PC, EasyOne Pro, and EasyOne Pro LAB |

Order Information

| Part number: | Product: |
|--------------|---|
| 2050-1 | Spirette box of 50 |
| 2050-5 | Spirette box of 200 |
| 2050-10 | Spirette box of 500 Not available in all countries |