

EasyOne Pro/LAB

Complete point-of-care pulmonary function testing

Making space for what matters in pulmonary health

The **EasyOne Pro/LAB** revolutionizes pulmonary function testing with its comprehensive range of tests, compact design, and user-friendly interface. Its capabilities, ranging from spirometry to lung volumes and DLCO, cover a wide spectrum of diagnostic needs in respiratory health.

The **EasyOne Pro/LAB** enables healthcare professionals to efficiently evaluate lung function, aiding in diagnosis, monitoring, and management of respiratory conditions. The combination of portability, ease of use, and accuracy makes the **EasyOne Pro/LAB** suitable for various clinical settings.



Precise performance



EasyOne Pro

- Spirometry
- DLCO
- Provocation testing

Lung volumes

- Single-breath helium dilution



EasyOne Pro LAB

- Spirometry
- DLCO
- Provocation testing

Lung volumes

- Multiple-breath nitrogen washout
- Single-breath helium dilution

Specifications

EasyOne Pro: FVC, FVL, SVC, MVV, Provocation and Pre/Post testing, DLCO

EasyOne Pro LAB: FVC, FVL, SVC, MVV, Provocation and Pre/Post testing, DLCO, MBW, LCI

Languages: Chinese, Croatian, Danish, Dutch, English, Finnish, French, German, Italian, Japanese, Norwegian, Polish, Portuguese, Russian, Spanish, Swedish, Turkish, Vietnamese

EMR connectivity: HL7, XML, GDT

Dimensions: 10.6 x 13.2 x 10.6", 17.65 lb
27 x 33.5 x 27 cm, 8 kg

Spirometry Standards: ATS/ERS 2022, 2019 & 2005, NIOSH, OSHA, SSA Disability

DLCO Standards: ATS/ERS 2017 & 2005

Article numbers

3000-1: EasyOne Pro

3100-1: EasyOne Pro LAB

2050-1: Spirette box of 50

2050-5: Spirette box of 200

2050-10: Spirette box of 500
Not available in all countries

3050-1: DLCO Barriettes box of 50

3050-2: DLCO Barriettes box of 100

3150-1: FRC Barriettes box of 40 (for use with EasyOne Pro LAB only)

3150-2: FRC Barriettes box of 80 (for use with EasyOne Pro LAB only)



Medical Technologies

n d d Medical Technologies
300 Brickstone Square, Suite 604 Andover, MA, 01810 USA
Phone: +1-978-470-0923

n d d Medizintechnik AG
Technoparkstrasse 1, CH-8005 Zurich, Switzerland
Phone: +41 44 512 65 00



Features

- Portable and lightweight design tailored for testing at the point of care.
- User-centric approach for intuitive operation, instant usability with no warm-up time.
- Hygienic solution with Spirette and Barriette consumables eliminates the risk of cross-contamination.
- Swiss precision engineering
- Minimal maintenance, allowing more time to be dedicated to patient care.

n d d TrueFlow

makes the difference

n d d's unique ultrasonic flow measurement is highly accurate in all flow ranges, independent of gas composition, pressure, temperature and humidity.

n d d **TrueFlow** is a resistance-free solution that does not require calibration during its lifetime.

n d d TrueCheck

automated precision

Patented method to verify performance of n d d machines

- Eliminates the need for gas-calibration as well as the need for an external DLCO simulator
- 5-point calibration including stability and linearity of tracer gas and CO sensors

EasyOne Connect

EasyOne Connect software offers one database and one platform for EasyOne point-of-care-solutions. Data management has never been easier.



ID: LP-0000048 Rev: 03
© Copyright 2025. All rights reserved