




Medical Technologies

Get to know Buddy,
your spirometry coach.



The future of spirometry

The **EasyOne Sky** transforms respiratory diagnosis and management with a fully connected spirometry solution.

Guided by Buddy, ndd's built-in spirometry coach, patients can perform tests more easily while providers achieve consistent, high-quality results. With secure cloud infrastructure protecting every test and real-time insights that give healthcare systems visibility into testing performance across their network, our ecosystem supports you from deployment to diagnosis, helping teams deliver better care with confidence.



EasyOne Sky

Guided spirometry. Connected insights.


EasyOne



Spirometry designed for your workflow.



Hi, I am Buddy,
your spirometry
coach!



- Buddy, your built-in spirometry coach, guides patients through each maneuver and provides real-time feedback to help achieve high-quality test results.
- Powered by nddCloud, **nddHub** is a web-based portal that provides centralized fleet and performance management across your entire network.
- Standardized setup, EMR integration, and structured support enable smooth spirometry rollout across any care setting.
- Training and continuous education with nddAcademy keep clinical teams confident and testing consistent over time.
- HITRUST-certified infrastructure automatically backs up and protects every test while enabling secure access anywhere.

Specifications

Test types: FVC, FVL, Pre/Post-testing

Languages: English, Spanish

Dimensions: 4.3" x 3.3" x 1.3", 4.9 oz / 111 x 84 x 33 mm, 140 g

Parameters: BEV, FEF₂₅, FEF₂₅₋₇₅, FEF₅₀, FEF₇₅, FET, FEV_{0.75}, FEV_{0.75}/FVC, FEV₁, FEV₁/FVC, FEV₁Q, FIVC, FVC, PEF, PIF, FVC %change, FEV1 %change

Predicted normal values spirometry: GLI (Quanjer) 2012, NHANES-III (Hankinson) 1999, Global GLI (Bowerman) 2023

Data management: EasyOne Mobile nddCloud nddHub

App compatibility: Windows

Article numbers

9500-1 EasyOne Sky Comfort subscription

9500-2 EasyOne Sky Enterprise subscription

9400-1 EasyOne Sky Comfort purchase

9400-2 EasyOne Sky Enterprise purchase

5050-50 EasyOne FlowTube box of 50

5050-200 EasyOne FlowTube box of 200



n d d Medical Technologies
300 Brickstone Square, Suite 604 Andover, MA, 01810 USA
Phone: +1 978 470 0923

n d d Medizintechnik AG
Technoparkstrasse 1, CH-8005 Zurich, Switzerland
Phone: +41 44 512 65 00



n d d TrueFlow
makes the difference

n d d's unique ultrasonic flow measurement is highly accurate in all flow ranges, independent of gas composition, pressure, temperature and humidity.

n d d **TrueFlow** is a hygienic and resistance-free solution that does not require calibration during its lifetime.

n d d Cloud
simplifying workflows

n d d's modern and secure connectivity solution is designed to connect to healthcare platforms and EHR/EMR systems. It is built on state-of-the-art cybersecurity standards.

SCAN ME

