

Patient Information

Name	Smith, Peter	Age	40 (11/8/1968)
ID	PSM-11213	Height	71.7 in
Gender	Male	Weight	176 lb
Ethnicity	Caucasian	BMI	24.1516711112268

Test Information

Test Type	FVC (ex only)	Predicted	NHANES III
Test Date	4/21/2009 6:43:49 PM	Value Selection	Best Value
Post Time		BTPS (IN/EX)	1.10/1.02
Interpretation	GOLD(2003)/Hardie		


This report was designed to meet the documentation requirements for Social Security Disability Determination.
 Calibration Data: (See Calibration Report Attached)

Smoker No

Test Result

Your FEV1 / Predicted: 91%

Parameter	Pred	LLN	Pre				%Pred
			Best	Trial 4	Trial 3	Trial 2	
FVC [L]	5.58	4.60	5.30	5.29	5.30	5.29	95
FEV1 [L]	4.43	3.60	4.02	4.02	3.96	3.91	91
FEV1/FVC	0.798	0.701	0.759	0.760	0.747	0.739	95
PEF [L/min]	633	487	709	701	703	709	112
BEV [L]	-	-	-	0.12	0.12	0.10	-
FET [s]	-	-	6.5	6.5	6.4	7.8	-
Session Quality	Pre	A (FEV1 Var=0.06L (1.6%); FVC Var=0.01L (0.1%))					



Volume = 10 mm per liter, Time = 20 mm per second