



Trust the ndd **EasyOne Connect** solution when it comes to your EMR integration project

Seamlessly integrate any **EasyOne** device into your EMR system using our **EasyOne Connect** software. **EasyOne Connect** software provides a fast, accurate and reliable connectivity solution to meet your needs.

✓ Advantages


- ✓ **Successful EMR integrations** with leading EMR vendors including Epic, Athena, Allscripts, eClinical Works, Cerner and more!
- ✓ **Robust user management solution** including active directory compatibility
- ✓ **Centralized patient database**
- ✓ **Support of all major data exchange protocols** including **HL7** and **API**
- ✓ **Dedicated connectivity team** for EMR implementation projects

Steps for a successful integration:



✓ It's that easy!

How does it work?



- ndd's EMR integration solution is **a file-based integration**.
- We utilize **SFTP or a network-based location** to communicate with an EMR vendor.
- The EMR vendor sends **EasyOne Connect** an **HL7 file** and then **EasyOne Connect** picks that file up and processes it. **That's it!**



Requirements from you (end customer)

All you need is a **suitable SFTP or network location** that ndd can utilize for communication with your EMR vendor.



Why the ndd solution?

- ✓ Hundreds of **successful integrations**
- ✓ We support your timeline, and can have an integration completed in **as little as one week**
- ✓ Our solution minimizes almost all sources of errors associated with data entry, to help you **deliver the best care to your patient**

Partnered with **over 30 different EMRs**

Ready to get started?

Submit a ticket to start your EMR project today!

