


Thermoprinter

Mobile Printer Solution
for EasyOne Air



TSC DA 200



TSC Alpha 4L

Connect your EasyOne Air via USB
to the thermoprinter.

Plug and Play - print without needing further installation.
With this solution you are even more mobile.

The reports print fast and are very easy to read.

The layout has been designed to be easily recognizable by other ndd products.

Model	TSC DA 200	TSC Alpha 4L
Resolution	8 dots/mm (203 DPI)	8 dots/mm (203 DPI)
Max. print speed	152.4 mm/second	102 mm/second
Max. print width	108 mm	104 mm
Dimension (WxHxD)	172 mm x 165 mm x 195 mm	160 mm x 191.6 mm x 79 mm
Weight	1.5 kg (3.3 lbs)	1.1 kg (2.43 lbs)
Interface	USB 2.0	USB 2.0
Power	AC 100-240V, 50-60Hz	7.2V DC, 5800 mAh Li-ion battery
Environmental conditions	5 ~ 40 °C, 10 ~ 85% non-cond.	-10 ~ 50 °C, 10 ~ 90% non-cond.

For further information and purchase conditions please contact TSC directly at www.tscprinters.com.
n d d distribution partners can purchase these thermoprinters under preferential conditions.