

Patient Information

Name PETER SMITH
ID NDD TECH
Age 50
Height 5 ft 4 in
Weight 190 lbs,BMI 32.8
Gender MALE
Ethnic CAUCASIAN
Smoker YES
Asthma NO

Test Information

Test Date/Time 03/20/2015 09:44am
Post Time --:--
Test Mode DIAGNOSTIC
Syst. Interpret. GOLD/Hardie
Predicted Ref Nhanes III
Value Select BEST VALUE
Tech ID NDD088
Automated QC ON
BTPS (IN/EX) ---/ 1.02

FVC Test Results Your FEV1 is 96% Predicted. Your Lung Age is 55

Parameter	Best	Pred	%Pred
FVC[L]	3.71	4.06	91
FEV1[L]	3.03	3.17	96
FEV1/FVC[%]	81.8	77.7	105
PEF[L/min]	573.9	509.2	113
FEF25-75[L/s]	3.09	2.92	106
FET[s]	8.72	--	--

Pre-Test FEV1 Var=0.05L 1.8%;FVC Var=0.06L 1.6%;Session Quality A
Syst. Interpret. Normal Spirometry

