

## EasyChecker Version History

This document lists all releases of the PC Software called EasyChecker. Some versions have been bundled with EasyWare and/or EasyWare Pro products. The changes though are cumulative, that means that changes introduced in a version are also included in later versions.

Version Date (dd.mm.yyyy)	Description
<b>V4.0.3.0 Rev. E</b> 15.07.2014	<b>Same EasyChecker Software, same EasyOne firmware archive but with new EasyChecker V4.0.3.0</b>
<b>V4.0.2.0 Rev. D</b> 14.06.2013	<b>Same EasyChecker Software, same EasyOne firmware archive but with new driver support for Windows 8 32-/64-Bit</b>
<b>V4.0.2.0 Rev. C</b> 31.01.2013	<b>Same EasyChecker Software but with new Firmware Files:</b> <ul style="list-style-type: none"> <li>• EZ1-STD-02250000-M-Res.zip replaces 0224</li> <li>• EZ1-STD-02250000-M-Res-Austria.zip replaces 0224</li> <li>• EZ1-STD-02250000-M-Res-Polish.zip replaces 0224</li> <li>• EZ1-STD-02250000-M-Res-Swedish-Norwegian.zip replaces 0224</li> <li>• EZ1-STD-02250000-P-Res.zip replaces 0224</li> <li>• EZ1-STD-02250000-P-Res-EE.zip replaces 0224</li> <li>• EZ1-STD-05100000-P-Res.zip replaces 0509</li> <li>• EZ1-STD-05100000-P-Res-EE.zip replaces 0509</li> <li>• EZ1-STD-06070000-P-Res.zip replaces 0606</li> <li>• EZ1-STD-06070000-P-Res-EE.zip replaces 0606</li> </ul>
<b>V4.0.2.0 Rev. B</b> 13.03.2012	<b>Same EasyChecker Software but with new Firmware Files:</b> <ul style="list-style-type: none"> <li>• EZ1-STD-05080000-P-Res.zip replaced by EZ1-STD-05090000-P-Res.zip</li> <li>• EZ1-STD-05080000-P-Res-EE.zip replaced by EZ1-STD-05090000-P-Res-EE.zip</li> <li>• EZ1-STD-06050000-P-Res.zip replaced by EZ1-STD-06060000-P-Res.zip</li> <li>• EZ1-STD-06050000-P-Res-EE.zip replaced by EZ1-STD-06060000-P-Res-EE.zip</li> </ul>
<b>V4.0.2.0</b> 21.06.2011  <b>Available on</b> <a href="http://www.ndd.ch">www.ndd.ch</a>  <b>“Download_01972_xx.exe”</b>	<b>First Version of n d d’s Device Software Checker which can be downloaded and installed independently from other n d d Software Products</b>  <b>Installer Software:</b> <ul style="list-style-type: none"> <li>• SmartInstaller.exe V3.0.2.0</li> </ul> <b>Bundled with:</b> <ul style="list-style-type: none"> <li>• DeviceCommDll.dll V1.6.9.0</li> <li>• Device Driver n d d USB Cradle for EasyOne (3.4.1.20) – WHQL approved – from Windows Update also</li> <li>• Device Driver n d d ScreenConnector for EasyOne (2.6) – WHQL approved – from Windows Update also</li> </ul> <b>Changes since previous version bundled with EasyWare V2.23.1:</b> <ul style="list-style-type: none"> <li>• Bug 50 Upgrade Problem around Austria 1994’ - under certain circumstances Austria 1994 Predicteds are erroneously classified as ‘must exist’.</li> </ul> <b>Firmware Files:</b> <ul style="list-style-type: none"> <li>• CRA-MRK-02050000.zip</li> <li>• CRA-STD-01040000.zip</li> <li>• CRD-STD-03031400.mot</li> <li>• EZ1-ALD-02150400-M-Res.zip</li> <li>• EZ1-ALD-02150400-P-Res.zip</li> <li>• EZ1-ANT-04020500-P-Res.zip</li> <li>• EZ1-CAA-04050000-P-Res.zip</li> <li>• EZ1-ESS-02220000-P-Res.zip</li> <li>• EZ1-GEC-04030000-P-Res-EU.zip</li> <li>• EZ1-GSK-02220000-P-Res.zip</li> <li>• EZ1-MCL-05050000-P-Res.zip</li> <li>• EZ1-MOT-02140400-M-Res.zip</li> <li>• EZ1-MOT-02140400-P-Res.zip</li> <li>• EZ1-MRK-03070000-P-Res.zip</li> <li>• EZ1-PF2-02170600-P-Res.zip</li> <li>• EZ1-SAV-02220000-P-Res.zip</li> <li>• EZ1-SIM-02220000-P-Res.zip</li> <li>• EZ1-STD-02220000-M-Res.zip</li> <li>• EZ1-STD-02220000-M-Res-Brasil.zip</li> <li>• EZ1-STD-02220000-M-Res-Polish.zip</li> <li>• EZ1-STD-02220000-M-Res-Swedish.zip</li> </ul>

# EasyWare Version History



File  
EasyCheckerVersionHistory-V05R.docx

Date  
15-July-14

Page  
2 of 2

<b>Version Date</b> (dd.mm.yyyy)	<b>Description</b>
	<ul style="list-style-type: none"><li>• EZ1-STD-02220000-P-Res-Polish.zip</li><li>• EZ1-STD-02220000-P-Res-Swedish-Norwegian.zip</li><li>• EZ1-STD-02240000-M-Res.zip</li><li>• EZ1-STD-02240000-M-Res-Austria.zip</li><li>• EZ1-STD-02240000-M-Res-Polish.zip</li><li>• EZ1-STD-02240000-M-Res-Swedish-Norwegian.zip</li><li>• EZ1-STD-02240000-P-Res.zip</li><li>• EZ1-STD-02240000-P-Res-EE.zip</li><li>• EZ1-STD-05050000-P-Res.zip</li><li>• EZ1-STD-05050000-P-Res-Polish.zip</li><li>• EZ1-STD-05080000-P-Res.zip</li><li>• EZ1-STD-05080000-P-Res-EE.zip</li><li>• EZ1-STD-06010000-P-Res.zip</li><li>• EZ1-STD-06050000-P-Res.zip</li><li>• EZ1-STD-06050000-P-Res-EE.zip</li></ul>