

Zürich, 17. March 2019

EasyWare Not Compatible with Win10 Notice

Dear EasyWare Customer,

The purpose of this notification is to inform you that EasyWare software is not compatible with the Windows 10 operating system. As all development activities on EasyWare have ceased, there will be no update to address this issue in the foreseeable future.

We would like to encourage you to upgrade your EasyOne device to the rechargeable, touch-screen operated [EasyOne Air](#). Otherwise you can also synchronize the EasyOne device with the [EasyOne Connect software](#) instead. If you encounter problems while exploring this option, please contact [ndd service and support](#).

Product Management
ndd Medizintechnik AG