

## EasyOne Mobile Software Version History

This document lists the most recent EasyOne Mobile Software releases. The changes made are cumulative, which means that changes introduced in one version are also included in later versions.

Version Date	Description
<b>V1.6.5</b> May 2026	Changes highlights since V1.5.6: <ul style="list-style-type: none"> <li>• Support for new development of nddHub:               <ul style="list-style-type: none"> <li>○ Multi-department support</li> <li>○ New role: Healthcare administrator</li> <li>○ Fleet management support</li> </ul> </li> <li>• Support for Cloud maintenance window</li> <li>• Extended patient profile with comment function and additional smoking information</li> <li>• Ability to edit first- &amp; last name of user profile</li> </ul>
<b>V1.5.6</b> August 2025	Changes highlights since V1.4.1: <ul style="list-style-type: none"> <li>• General Changes               <ul style="list-style-type: none"> <li>○ Separate workflow for Multi-flow CalCheck</li> <li>○ Add trial to already closed session within 24h</li> <li>○ Show PDF report before saving</li> <li>○ Updating historic test results after input parameter changed</li> <li>○ Enabling push-notifications</li> <li>○ Synchronization improvements</li> <li>○ Introduction of Blazor Hybrid</li> <li>○ Delete Order functionality</li> </ul> </li> <li>• Bug Fixes               <ul style="list-style-type: none"> <li>○ Added missing time zone information</li> </ul> </li> </ul>
<b>V1.4.1</b> January 2025	Changes highlights since V1.3.4: <ul style="list-style-type: none"> <li>• General Changes               <ul style="list-style-type: none"> <li>○ Sharing patients between users within the same organization</li> <li>○ Change of user interface when entering patients date of birth</li> <li>○ New default settings upon initial setup:                   <ul style="list-style-type: none"> <li>▪ Show incentive</li> <li>▪ GLI 2021</li> </ul> </li> </ul> </li> <li>• Bug Fixes               <ul style="list-style-type: none"> <li>○ Made e-mail comparison during sign-in/sign-up case-insensitive</li> <li>○ Fixed issue with “Forgot your password?”-Workflow</li> <li>○ Fixed wrong height &amp; weight conversion when changing language from English to Spanish</li> </ul> </li> </ul>
<b>V1.3.4</b> November 2024	First official distributed release.